



856 Antares

Cervical collar set 6 pieces with carrying case



FEATURES

ANTARES cervical collar is an highly professional instrument indispensable for the emergency operator. It is composed by 2 pieces to be put on the front and back part of the neck enabling thus the perfect immobilization of the rachis. Central hole for tracheotomy and easy positioning thanks to the rip-off fastening which different colours indicate the various sizes. It is made with transparent polyethylene with E.V.A. washable inner lining fixed by plastic rivets. Available in 6 different sizes.

MAIN FEATURES

Usage	Immobilization
Trauma	

SEE ALSO



Sirius 874

Cervical collar set 4 pieces with carrying case Art. 862 - 864 - 866 - 868



Sirius 876

Cervical collar set 6 pieces with carrying case



Sirius 874

Cervical collar set 4 pieces with carrying case Art. 862 - 864 - 866 - 868



Antares 854

Cervical collar set 4 pieces with carrying case Art. 842 - 844 - 846 - 848



Antares 840

Rigid cervical collar in 2 pieces size INFANT



Antares 842

Rigid cervical collar in 2 pieces size PEDIATRIC



Antares 844

Rigid cervical collar in 2 pieces size STOUT



Antares 846

Rigid cervical collar in 2 pieces size STANDARD



Antares 848

Rigid cervical collar in 2 pieces size MEDIUM



Antares 850

Rigid cervical collar in 2 pieces size TALL