



874 Sirius

Cervical collar set 4 pieces with carrying case
Art. 862 - 864 - 866 - 868



FEATURES

Rigid cervical collar in 1 piece made with transparent polyethylene with E.V.A. washable inner lining fixed by plastic rivets. Adjustable rip-off fastening which different colours indicate the various sizes. Central hole for tracheotomy. Available in 6 different sizes.
Art. 862 - 864 - 866 - 868

MAIN FEATURES

Usage

Trauma

Immobilization



SEE ALSO



Sirius 876
Cervical collar set 6 pieces with
carrying case



Antares 856
Cervical collar set 6 pieces with
carrying case



Antares 854
Cervical collar set 4 pieces with
carrying case Art. 842 - 844 - 846 -
848



Sirius 862
Rigid cervical collar in 1 piece size
PEDIATRIC



Sirius 864
Rigid cervical collar in 1 piece size
STOUT



Sirius 866
Rigid cervical collar in 1 piece size
STANDARD



Sirius 868
Rigid cervical collar in 1 piece size
MEDIUM