



884 Halley

Vacuum splint for forearm



FEATURES

Produced with absolute reliability materials this splints allow the perfect immobilization of the lower and upper limbs by taking away the inside air. They are supplied with suction valve positioned on the out-side to avoid problems to the patient. It is possible to order the single vacuum splint or the complete set including splints for leg, arm, forearm, little pump for the vacuum and nylon bag. All splints are X-ray transparent.

MAIN FEATURES

Size	Adults
Usage	Immobilization
Trauma	
Width (cm)	34
Material	PVC
Length (cm)	52

IMAGES



SEE ALSO



Antares 854

Cervical collar set 4 pieces with carrying case Art. 842 - 844 - 846 - 848



Halley 882

Vacuum splint for arm



Halley 880

Vacuum splint for leg



Halley 886

Vacuum splint set with pump and carrying case

ACCESSORIES



889

Little aluminium pump
