



9012 Galaxy

"GALAXY" Orange spinal board



FEATURES

Spinal board made of plastic material at high strength. Designed paying particular attention to the difficulties met during every day's interventions; they have 4 holes for the quick and total fixing of the head immobilizer and two cavities where the board lays on the floor, when the base is blocked in the traditional way, that allow to avoid damages to the rip-off straps during the usage of the spinal board or the accommodation in the ambulance. It has 12 handles for the transport and is supplied with 3 belts with rapid unhooking buckle (Art. 605/MEB). X-ray translucent and available in 3 different colours: yellow, orange and blue.

Device certified according to the European harmonized safety standards UNI EN 1865.

MAIN FEATURES

Туре	
Size	Adults
Usage	Immobilization
Length (cm)	183
Width (cm)	41
Thickness (cm)	5
Weight (kg)	5,5
Loading capacity (kg)	150
Trauma	₩
Material	ABS

SEE ALSO



ACCESSORIES



Spider 690Belts immobilization system Spider for spinal board



Jelly 629 Red standard head immobilizer for spinal board



12144/N Carrying case for spinal board

STANDARD ACCESSORIES



605/MEB50 mm orange belt in one piece with plastic buckle and ME.BER. ribbon