



## 910 Winner

---

Self loading stretcher



---

### FEATURES

Article 910 has been designed to offer a practical and safe instrument to the operator in case of emergency. The stretcher, built with anodized aluminium, is strong and handy thanks to its four wheels diam. 150, two fixed and with brakes and two rotating. Important feature of this article is the loading height regulation system. The automatic loading stretcher "WINNER" in fact is supplied with two sturdy perforated plates that allow to change the little front wheels height in order to load the stretcher on ambulances of various model and with different inside fittings.

- Back rest and lying part made with thermoplastic material reinforced to allow cardiac massage. Easy to clean and to disinfect
- Adjustable back rest 0° - 80° moved by a gas piston which allow to fix the back rest safely in any position
- Side protections completely overturnable
- Blocking device for folded legs with a system that avoid the compression of the blocking bar
- Safety lever for the legs positioned near the unlocking device allowing thus the release operation for the loading, keeping the hands on the stretcher
- Vertical legs protected by nylon wedges
- Automatic centring device mounted on rotating wheels. This system automatically block the back wheels in the central position during the loading of the stretcher on the ambulance without having to turn the wheels manually

---

## MAIN FEATURES

---

Size	Adults
Usage	Transport
Lying surface	Without Trendelenburg
Wheel type	2 Rotating Wheels + 2 Small Fixed Wheels
Type	Self loading stretcher
Height (cm)	87,5
Width (cm)	56
Weight (kg)	34
Loading capacity (kg)	170
Folded height (cm)	26
Max loading height (cm)	65
Length (cm)	194
Min loading height (cm)	45
Material	Anodized aluminium
Length with folded trolley (cm)	162

---

## SEE ALSO

---



**Winner 910 PROOF**  
Self loading stretcher - certified

---

## ACCESSORIES



**Front blocking device**  
Front blocking device for automatic loading stretcher



**Rear fixing system**  
Rear fixing system for roll-in H=148 mm



**950/2**  
Rear fixing system for roll-in H=160 mm



**950/1**  
Rear fixing system for roll-in H=188 mm



**936**  
Phleboclysis stand for automatic loading stretcher Winner and Frog



**939**  
Instrument support for stretchers



**938**  
Oximeb oxygen cylinder holder



**612/AN-2/MEB**  
Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



**612/N-2/MEB**  
Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



**608/A-MEB**  
2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



**608/N-MEB**  
2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



**12025**  
Black multi usage object holder for automatic loading stretcher



**7013/N**  
3 pcs. Skay fire-proof electrowelded mattress



**7014/A**  
4 pcs. Skay fire-proof electrowelded mattress



**12018N**  
Telescopic foldable phleboclysis stand for self loading stretchers



**12070/BG**  
"COVER OX 2-3 BG"- professional cylinder holder for 2 and 3 lt oxygen bottles



**12070/BR**  
COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



**12054**  
Oximeb Maxi - Oxygen cylinder holder



**12072/BC**  
"COVER OX 5 BC"- professional  
cylinder holder for 5 lt oxygen  
bottles



**12072/BR**  
"COVER OX 5 BR"- professional  
cylinder holder for 5 lt oxygen  
bottles

---