



9100 Bunny

Pediatric vacuum mattress with independent chambers



FEATURES

Mattress made of polyurethane with fire retardant treatment, 14 independent chambers for a perfect and safe immobilization. Usable between -30°C and +70°C. Extreme light weight and strong, 8 ergonomic handles padded and reinforced, safety belts and a perfect containment of the patient, give this mattress high performances. Usable with X-ray, CR and MRI. Weight: 3 kg.

MAIN FEATURES

Size	Children
Usage	Immobilization
Length (cm)	120
Width (cm)	70
Weight (kg)	3
Pediatric	<input checked="" type="checkbox"/>
Trauma	<input checked="" type="checkbox"/>
Material	Polyurethane

SEE ALSO



Bunny Plus 9102

Independent chambers vacuum mattress pediatric with integrated head immobilizer

ACCESSORIES



888
Aluminium aspiration pump



896/DF
Abs aspiration pump



9100-01
Carrying bag for pediatric vacuum mattress

STANDARD ACCESSORIES



894-001
Aspiration valve for vacuum splints
and mattress
