



9102 Bunny Plus

Indipendent chambers vacuum mattress pediatric with integrated head immobilizer



FEATURES

Mattress with integrated head immobilizer to allow a fast head immobilization. Made of polyurethane treated with fire retardant, **14 independent chambers** for a perfect and safe immobilization. Usable between - 30°C and + 70°C. Extreme light weight and strenght, 8 ergonomic handles padded and reinforced, safety belts and a perfect containment of the patient, give this mattress high performances. Usable with X-ray, CT and MRI.

MAIN FEATURES

Size	Children
Usage	Immobilization
Length (cm)	120
Width (cm)	70
Weight (kg)	3,5
Pediatric	<input checked="" type="checkbox"/>
Trauma	<input checked="" type="checkbox"/>
Material	Polyurethane

SEE ALSO



Bunny 9100
Pediatric vacuum mattress with
independent chambers

ACCESSORIES



888
Aluminium aspiration pump



896/DF
Abs aspiration pump

STANDARD ACCESSORIES



894-001
Aspiration valve for vacuum splints
and mattress
