



936

Phleboclysis stand for automatic loading stretcher Winner and Frog



FEATURES

Telescopic phleboclysis stand for stretcher. Thanks to its fixing system it is possible to fix it in every point on the stretcher's frame. **Dimensions: minimum height: cm 47,5 maximum height: cm 130,5**

MAIN FEATURES

Min height (cm)	47,5
Max height (cm)	130,5
Weight (kg)	0,5
Material	Aluminium alloy

IMAGES



SEE ALSO



936/A

Phlebotoclysis stand for stretchers
Art. 1510 Proof and Art. 1515 Proof



12018N

Telescopic foldable phlebotoclysis
stand for self loading stretchers

ACCESSORY OF



7265 Proof

Aluminium self loading stretcher
"Frog lite Advance" with
exctractable handles and
Trendelenburg