



## 936/A

Phleboclysis stand for stretchers Art. 1510 Proof and Art. 1515 Proof



### FEATURES

Telescopic phleboclysis stand for stretcher. Thanks to its fixing system it is possible to fix it in every point on the stretcher's frame. **Dimensions: minimum height: cm 47,5 maximum height: cm 130,5**

### MAIN FEATURES

Min height (cm)	47,5
Max height (cm)	130,5
Weight (kg)	0,5
Material	Aluminium alloy

### IMAGES



---

SEE ALSO

---



**936**

Phlebotomy stand for automatic loading stretcher Winner and Frog



**12018N**

Telescopic foldable phlebotomy stand for self loading stretchers

---