



## 965 Front blocking device

Front blocking device for automatic loading stretcher



### FEATURES

Blocking device to be fixed on the ambulance platform that allows to guide and block the automatic loading stretcher in the front part. Made of stainless steel.

### MAIN FEATURES

|             |      |
|-------------|------|
| Length (cm) | 51,5 |
| Width (cm)  | 48   |
| Height (cm) | 10   |
| Weight (kg) | 1,4  |

### SEE ALSO



**Winner 914**  
Self loading stretcher with  
trendelenburg and fowler



**Winner 910**  
Self loading stretcher