



7122/AZ Gecko

“GECKO” self loading stretcher entry-level light blue



FEATURES

Studied to optimize quality-price ratio. Build with painted stainless steel. This device grants all features of other self loading stretchers as: adjustable back rest, 2 fixed wheels with brake and 2 rotating wheels Ø 200 mm self centring, folding side protections. Supplied with belts.

MAIN FEATURES

| | |
|-------------------------|-----------------------|
| Type | |
| Usage | Transport |
| Lying surface | Without Trendelenburg |
| Height (cm) | 90,5 |
| Max loading height (cm) | 62 |
| Min loading height (cm) | 62 |
| Folded height (cm) | 30 |
| Width (cm) | 56 |
| Length (cm) | 197 |
| Weight (kg) | 37 |
| Loading capacity (kg) | 170 |
| Material | Painted steel |

IMAGES



SEE ALSO



Gecko 7122/V
"GECKO" self loading stretcher
entry-level green

ACCESSORIES



939
Instrument support for stretchers



Rear fixing system
Rear fixing system for roll-in H=148 mm



Front blocking device
Front blocking device for automatic loading stretcher



938
Oximeb oxygen cylinder holder



7016/A
PVC orange electrowelded mattress (one-piece)



612/AN-2/MEB
Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



612/N-2/MEB
Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



12080-001
Monting kit for bariatric platform on self loading stretcher



12010
Phleboclysis stand for stretchers Mercury and Gecko



12025
Black multi usage object holder for automatic loading stretcher



608/A-MEB
2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



608/N-MEB
2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



12080
"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher

STANDARD ACCESSORIES



603/MEB
Nylon bag containing 3 belts Art. 607/Ameb