



# 7122/V Gecko

“GECKO” self loading stretcher entry-level green



## FEATURES

Studied to optimize quality-price ratio. Build with painted stainless steel. This device grants all features of other self loading stretchers as: adjustable back rest, 2 fixed wheels with brake and 2 rotating wheels Ø 200 mm self centring, folding side protections. Suplied with belts.

## MAIN FEATURES

Type	Transport
Usage	Without Trendelenburg
Lying surface	
Height (cm)	90,5
Max loading height (cm)	62
Min loading height (cm)	62
Folded height (cm)	30
Width (cm)	56
Length (cm)	197
Weight (kg)	37
Loading capacity (kg)	170
Material	Painted steel

---

IMAGES

---



---

SEE ALSO

---



**Gecko 7122/AZ**  
"GECKO" self loading stretcher  
entry-level light blue

---

## ACCESSORIES



**7016/A**  
PVC orange electrowelded  
mattress (one-piece)



**939**  
Instrument support for stretchers



**Rear fixing system**  
Rear fixing system for roll-in H=148  
mm



**Front blocking device**  
Front blocking device for  
automatic loading stretcher



**938**  
Oximeb oxygen cylinder holder



**612/AN-2/MEB**  
Adjustable thoracic abdominal belt  
orange/black with MeBer band and  
4 points fastening for self loading  
stretchers and sedan chairs



**612/N-2/MEB**  
Adjustable thoracic abdominal belt  
black with MeBer band with 4  
points fastening for self loading  
stretchers and sedan chairs



**12080-001**  
Monting kit for bariatric platform  
on self loading stretcher



**12025**  
Black multi usage object holder for  
automatic loading stretcher



**12080**  
"MPB" wide bariatric platform for  
self-loading stretchers Frog  
Mercury and Gecko complete of  
belts with reel and fixing for  
stretcher



**12010**  
Phleboclysis stand for stretchers  
Mercury and Gecko



**608/A-MEB**  
2 pieces orange belt with metallic  
hook and MeBer band (L. 155 cm)



**608/N-MEB**  
2 pieces black belt with metallic  
hook and MeBer band (L. 155 cm)



**12070/BG**  
"COVER OX 2-3 BG"- professional  
cylinder holder for 2 and 3 lt oxygen  
bottles



**12070/BR**  
COVER OX2-3 BR professional  
cylinder holder with belt and  
orange MeBer logo for 2 and 3 lt  
oxygen bottles



**12072/BG**  
"COVER OX 5 BG"- professional  
cylinder holder for 5 lt oxygen  
bottles



**12072/BR**

"COVER OX 5 BR"- professional  
cylinder holder for 5 lt oxygen  
bottles

---

---

## STANDARD ACCESSORIES

---



**603/MEB**

Nylon bag containing 3 belts Art.  
607/Ameb

---