



431

Arcos rigid splint set with carrying bag



FEATURES

New rigid splints, made with fire-proof and water proof PVC, allow a correct immobilization of the limbs. They are X-ray translucent. Equipped with 4 Velcro® closures with adjustable length allowing operator to apply them easily and safely. Shaped to strongly and easily fix to ankle and upper part of the leg, they are supplied with foot support integrated.

Available in two sizes:

Adult arm/ pediatric leg: length cm. 52 - base cm 9

Adult leg: length cm 74 - base cm 13

Pediatric limb: length 30 cm, base 9 cm

MAIN FEATURES

Trauma

Usage



Immobilization

IMAGES



SEE ALSO



425

Arcos long rigid splint for lower limbs



427

Arcos short rigid splint for upper limbs



429

Arcos Splint for pediatric limbs
